


Endoscopic Needle

For use with Endoscopic Marker

The success of an endoscopic tattoo highly depends on the application of the endoscopic marker at the correct depth. This needle is designed for safe injection of the endoscopic marker at the correct depth.

Features:

- Highly lubricant PTFE outertube for easy passage through the working channel.
- Strong outertube facilitates easy advancement through complicated anatomic structures.
- Very sharp needle tip specifically designed for safe and easy submucosal penetration.
- Shielded tip prevents the needle from damaging the outertube or the working channel.



Artikel nr.	Omschrijving	Aantal
EMN2504	Endoscopic Marker Injection Needle	10 stuks